


① 肌氨酸 (M.W.89.09)  
Sarcosine



② 丙氨酸 (M.W.89.09)  
Alanine

### 色谱条件

色谱柱	: CAPCELL PAK CR 1:4 S5; 2.0 mm i.d.×150 mm
流动相	: 100mM 甲酸铵, 0.1%甲酸/乙腈=15/85
流速	: 0.2 mL/min
温度	: 40 °C
检测	: MS ESI Positive, MRM 90.1/44.1Da
进样量	: 1 μL
样品浓度	: 1ppm, 溶剂为流动相